

INSTRUCTIONS FOR USE

VITROS Immunodiagnostic Products SARS-CoV-2 Antigen Extraction Buffer

REF 619 9944

Rx ONLY

Intended Use

For *in vitro* diagnostic and laboratory professional use.

For use to extract the viral load from swab samples placed in appropriate transport media for use on the VITROS 3600 Immunodiagnostic System and the VITROS 5600/XT 7600 Integrated Systems.

Warnings and Precautions

WARNING: *Potentially Infectious Material*

Handle, use, store and dispose of solid and liquid waste from samples and test components, in accordance with procedures defined by appropriate national biohazard safety guideline or regulation (e.g. CLSI document M29).¹

WARNING: *Contains Mixture, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone (CAS 55965-84-9)²*

The VITROS SARS-CoV-2 Antigen Extraction Buffer contains ≥ 0.0015 – $< 0.06\%$ of Mixture, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone. H317: May cause an allergic skin reaction. P280: Wear protective gloves. P333 + P313: If skin irritation or rash occurs: Get medical advice/attention. P362 + P364: Take off contaminated clothing and wash before reuse.

Refer to www.orthoclinicaldiagnostics.com for the Safety Data Sheets and for Ortho contact information.

WARNING



Materials Provided

1 extraction buffer pack containing:

- 4 bottles of VITROS SARS-CoV-2 Antigen Extraction Buffer (28 mL) with antimicrobial agent

Materials Required but Not Provided

- Appropriate volume pipette and sample containers for extraction

Refer to the VITROS Immunodiagnostic Products SARS-CoV-2 Antigen Reagent Pack and Calibrator instructions for use.

Storage, Preparation and Handling

- The VITROS Immunodiagnostic Products SARS-CoV-2 Antigen Extraction Buffer is supplied ready to use.
- Store unopened at 2–8 °C (36–46 °F).

INSTRUCTIONS FOR USE

Testing Procedure

- Do not freeze.
- Use opened vial within 4 weeks of opening.
- Do not use beyond the expiration date.

Testing Procedure

Refer to the VITROS Immunodiagnostic Products SARS-CoV-2 Antigen Reagent Pack and Calibrator instructions for use.

Note: Do not use visibly damaged product.




































Not all products and systems are available in all countries.

References

1. CLSI. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline - Fourth Edition*. CLSI document M29-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
2. Regulation (EC) No 1272/2008 of the European Parliament and of the Council of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

Glossary of Symbols

The following symbols may have been used in the labeling of this product.

	Do Not Reuse		Upper Limit of Temperature		Range
	Use by or Expiration Date (Year-Month-Day)		Lower Limit of Temperature		Range of Means
	Batch Code or Lot Number		Temperature Limitation		Midpoint
	Serial Number		Consult Instructions for Use		Revised
	Catalog Number or Product Code		Attention: The Instructions for Use (IFU) has been updated		Supersedes
	Caution		For use in Slide Supply 1		Contains Sufficient for "n" Tests
	Keep Dry (Protect from Moisture/Humidity)		For use in Slide Supply 2		<i>in vitro</i> Diagnostic Medical Device
	Manufacturer		SI Units		Der Grüne Punkt (the Green Dot). Manufacturer follows certain packaging material waste disposal management regulations
	Date of Manufacture		Conventional Units		Estimated within-lab SD
	Authorized Representative in the European Community		Value		Serious Health Hazards
	Corrosive		Flammable		Environmental or Aquatic Toxicity
	Health Hazards		Acute Toxicity		

INSTRUCTIONS FOR USE

Revision History

Revision History

Date of Revision	Version	Description of Technical Changes*
2020-10-22	1.0	Initial version of document

* The change bars indicate the position of a technical amendment to the text with respect to the previous version of the document.


When this Instructions For Use is replaced, sign and date below and retain as specified by local regulations or laboratory policies, as appropriate.

Signature

Obsolete Date

Conditions of supply: all supplies are made subject to the standard terms and conditions of Ortho Clinical Diagnostics or its distributors. Copies of these are available on request.

Distributed in the US by:
Ortho-Clinical Diagnostics, Inc.
100 Indigo Creek Drive
Rochester, NY 14626



Ortho-Clinical Diagnostics
1500 Boulevard Sébastien Brant
B.P. 30335
67411 Illkirch
CEDEX, France



Ortho-Clinical Diagnostics
Felindre Meadows
Pencoed
Bridgend
CF35 5PZ
United Kingdom

VITROS is a trademark of Ortho Clinical Diagnostics.
© Ortho Clinical Diagnostics, 2020

Ortho Clinical Diagnostics