

INSTRUCTIONS FOR USE

hsTnI

VITROS Immunodiagnostic Products hs Troponin I Reagent Pack

REF 684 4436

VITROS Immunodiagnostic Products hs Troponin I Calibrators

REF 684 4437

Intended Use

For *in vitro* diagnostic use only.

VITROS Immunodiagnostic Products hs Troponin I Reagent Pack

For the quantitative measurement of cardiac troponin I (cTnI) in human serum and plasma (heparin) using the VITROS ECI/ECiQ/3600 Immunodiagnostic Systems and the VITROS 5600/XT 7600 Integrated Systems.

Cardiac troponin I elevations are associated with myocardial damage, and in patients with acute coronary syndrome (ACS), troponin values are used to aid in the diagnosis of myocardial infarction (MI). The test is further indicated for risk stratification of mortality, myocardial infarction or coronary revascularization in patients with acute coronary syndrome.

VITROS Immunodiagnostic Products hs Troponin I Calibrators

For use in the calibration of the VITROS ECI/ECiQ/3600 Immunodiagnostic Systems and the VITROS 5600/XT 7600 Integrated Systems for the quantitative measurement of cardiac troponin I in human serum and plasma (heparin).

Summary and Explanation of the Test

Troponin I, troponin T, and troponin C are components of the contractile apparatus in cardiac and skeletal muscle. Troponin I (TnI) in conjunction with troponin T and troponin C, regulates the calcium dependent interaction of actin and myosin.¹

Specificity of cTnI for myocardium: Two of the three isoforms of TnI identified are expressed in the skeletal muscle, while the third one is selectively expressed in the cardiac muscle (myocardium).^{2,3} This cardiac TnI isoform (cTnI) has an additional 31 amino acid residues at the N terminus, not found in skeletal troponin.⁴ Consequently, increases in cTnI values have not been reported to occur following injury to non-cardiac tissues and are considered specific for myocardial damage.^{5,6} Other biomarkers, e.g. creatine kinase MB isoform (CKMB), are less sensitive and less specific.⁷ The VITROS hs Troponin I (hsTnI) test uses monoclonal antibodies specifically directed against human cTnI.

cTnI and Myocardial Injury: Elevated cardiac troponin (cTn) values above the 99th percentile Upper Limit of Normal (URL) indicate the presence of myocardial injury, which can be acute, with dynamic changes in cTn levels, or chronic, with persistently elevated cTn levels.⁵ Conditions resulting in acute or chronic myocardial injury and elevated cTnI levels include but are not limited to heart failure, myocarditis, cardiomyopathy, critical illness, cardiotoxic chemotherapy, pulmonary embolism, pulmonary hypertension, stroke, sepsis and chronic kidney disease.⁵ These cardiac and systemic conditions should also be considered when interpreting results from any cTnI test method.

cTnI values provide important prognostic information about the short and long term risk of cardiovascular death.⁸⁻¹¹ The assessment of the prognosis can be useful in identifying patients most likely to benefit from specific therapeutic interventions.

cTnI and Myocardial Infarction (MI): An acute MI is defined by an acute myocardial injury associated with a rising and/or falling pattern of cTn values with at least one value above the 99th percentile URL in the clinical context of myocardial ischaemia (i.e. characteristic symptoms, ECG changes, imaging evidence).⁵ Clinical studies have demonstrated that cTnI is detectable in the bloodstream 4-6 hours after an acute MI and remains elevated for several days thereafter.¹²

Importance of cTnI Serial Measurement: Guidelines state that serial cTnI measurement is needed to differentiate an acute MI from other conditions resulting in cTnI elevations.^{5,13} cTnI release depends on blood flow, therefore, patient to patient variability is expected in the time to peak value, the time when cTnI exceeds the 99th percentile, or when changing pattern of values can be observed. Around peak levels (late presenters) and on the downslope curve (very late presenters), the change in cTnI values might not be as significant as on the upslope (early presenters).⁵

Algorithms incorporating cTnI to Rule-in or Rule-out MI vary by region and clinical guidelines. However, they should be used in conjunction with all available clinical information (i.e. detailed assessment of chest pain, ECG, imaging). In patients presenting very early a second cardiac troponin level should be obtained later, due to the time dependency of troponin release; also, additional serial cardiac troponin testing should be pursued if the clinical suspicion remains high or whenever the patient develops recurrent chest pain.¹³

Principles of the Procedure


An immunometric immunoassay technique is used. Cardiac troponin I present in the sample reacts simultaneously with streptavidin-conjugated antibody (mouse monoclonal anti-cTnI), bound by biotin-BSA on the wells, and a horseradish peroxidase (HRP)-labeled antibody conjugate (mouse monoclonal anti-cTnI). The antigen-antibody complex is captured by the antibody on the wells. Unbound materials are removed by washing. The bound HRP conjugate is measured by a luminescent reaction.¹⁴ A reagent containing luminogenic substrates (a luminol derivative and a peracid salt) and an electron transfer agent is added to the wells. The HRP in the bound conjugate catalyzes the oxidation of the luminol derivative, producing light. The electron transfer agent (a substituted acetanilide) increases the level of light produced and prolongs its emission. The light signals are read by the system. The amount of HRP conjugate bound is directly proportional to the concentration of cardiac troponin I present.

Test Type	System *	Incubation Time	Time to first result	Test Temperature	Reaction Sample Volume
Immunometric	ECi/ECiQ, 3600, 5600, XT 7600	8 minutes	15 minutes	37 °C	60 µL

* Not all products and systems are available in all countries.

Reaction Scheme

Biotin-BSA
Coated Well



Warnings and Precautions

WARNING:

Potentially Infectious Material

Human blood products provided as components of the VITROS hs Troponin I Calibrators have been obtained from donors who were tested individually and who were found to be negative for hepatitis B surface antigen, and for antibodies to human immunodeficiency virus (HIV 1+2) and hepatitis C virus (HCV), using approved methods (enzyme immunoassay). Treat as if capable of transmitting infection.

Use caution when handling material of human origin. Consider all samples potentially infectious. No test method can offer complete assurance that hepatitis B virus, HCV, HIV 1+2 or other infectious agents are absent. Handle, use, store and dispose of solid and liquid waste from samples and test components, in accordance with procedures defined by appropriate national biohazard safety guideline or regulation (e.g. CLSI document M29).¹⁵

WARNING:

Contains ProClin 950 (CAS 2682-20-4)

The VITROS hs Troponin I assay reagent contains 1.0% ProClin 950. The VITROS hs Troponin I HRP conjugate reagent contains 0.5% ProClin 950. H315: Causes skin irritation. H317: May cause an allergic skin reaction. H319: Causes serious eye irritation. P280: Wear protective gloves/protective clothing/eye protection/face protection. P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P337 + P313: If eye irritation persists: Get medical advice/attention. P302 + P352: IF ON SKIN: Wash with plenty of soap and water. P333 + P313: If skin

irritation or rash occurs: Get medical advice/attention. P363: Wash contaminated clothing before reuse.

Refer to www.orthoclinical.com for the Safety Data Sheets and for Ortho contact information.

WARNING



Reagents

Reagent Pack Contents

- 1 reagent pack containing:
 - 100 coated wells (biotin-BSA; streptavidin-mouse monoclonal anti-cTnI, 4 µg/mL).
 - 8.2 mL assay reagent (buffer with horse serum, bovine gamma globulin, bovine serum albumin, and antimicrobial agent).
 - 7.0 mL conjugate reagent (HRP-mouse monoclonal anti-cTnI, 5 µg/mL, in buffer with bovine serum albumin and antimicrobial agent).

Reagent Pack Handling

- The reagent pack is supplied ready for use.
- The reagent pack contains homogeneous liquid reagents that do not require shaking or mixing prior to loading on the system.
- As with all immunoassay protein-based solutions, inappropriate handling of the reagent pack can cause foam to occur on the surface of the reagent. Avoid agitation, which may cause foaming or the formation of bubbles.
 - If reagent packs are dropped or agitated, small levels of fine foam could be generated that may not be detected by the system.
 - Reagent packs containing fine foam that is not detected by the system, may show a negative bias.
- If you must use a dropped or agitated reagent pack before it has been allowed to settle, you should verify performance by running high and low quality control samples in duplicate after loading the pack on the system.

Reagent Pack Storage and Preparation

Reagent	Storage Condition		Stability
Unopened	Refrigerated	2–8 °C (36–46 °F)	expiration date
Opened	On system	System turned on	≤4 weeks
Opened	Refrigerated	2–8 °C (36–46 °F)	≤4 weeks

- The VITROS hs Troponin I Reagent Pack is suitable for use until the expiration date on the carton when stored and handled as specified. Do not use beyond the expiration date.
- Do not freeze unopened reagent packs.
- Load reagent packs directly from refrigerated storage to minimize condensation.
- Opened reagent packs are moisture/humidity sensitive. Store opened refrigerated reagent packs in a sealed VITROS Immunodiagnostic Products Reagent Pack Storage Box with desiccant.

Calibrator Contents

- 1 set of VITROS hs Troponin I Calibrators 1, 2 and 3 (2 mL, human cTnI in human serum/buffer matrix containing EDTA with antimicrobial agents); nominal values 0; 40 and 27,000 ng/L (pg/mL)
- Lot calibration card
- Protocol card
- 24 calibrator bar code labels (8 for each calibrator)

Calibrator Handling

- Use only with reagent packs of the same lot number. Mix thoroughly by inversion and bring to 15–30 °C (59–86 °F) before use. Each pack contains sufficient volume for a minimum of 8 determinations of each calibrator.

- Handle calibrators in stoppered containers to avoid contamination and evaporation. To avoid evaporation, limit the amount of time calibrators are on the system. Refer to the operating instructions for your system. Return to 2–8 °C (36–46 °F) as soon as possible after use, or load only sufficient volume for a single determination.

Calibrator Storage and Preparation

Calibrator	Storage Condition	Stability
Unopened	Frozen ≤-20 °C (-4 °F)	expiration date
Opened	Refrigerated 2–8 °C (36–46 °F)	≤2 Days
Opened	Frozen ≤-20 °C (-4 °F)	≤13 Weeks

- VITROS hs Troponin I Calibrators are supplied frozen.
- VITROS hs Troponin I Calibrators are suitable for use until the expiration date on the carton when stored and handled as specified. Do not use beyond the expiration date.
- Allow vials to stand upright at room temperature, for a minimum of 30 minutes to thaw. Thoroughly mix calibrators by inversion before use. Do not shake.
- Once thawed and opened calibrators may be stored at 2–8 °C (36–46 °F) for up to 2 days or frozen (with no more than 1 freeze-thaw cycle) for up to 13 weeks.
- The VITROS hs Troponin I test uses 60 µL of calibrator for each determination. Transfer an aliquot of each calibrator into a sample container (taking account of the minimum fill volume of the container), which may be bar coded with the labels provided. For details on minimum fill volume of sample cups or containers, refer to the operating instructions for your system.

Specimen Collection, Preparation and Storage

Patient Preparation

No special patient preparation is necessary.

Specimens Recommended

- Serum
- Plasma (Lithium Heparin)

Specimens Not Recommended

- Do not use turbid specimens. Turbidity in specimens may affect test results.
- Do not use hemolysed samples as hemolysis may affect test results.

Special Precautions

IMPORTANT:

Certain collection devices have been reported to affect other analytes and tests.¹⁶ Owing to the variety of specimen collection devices available, Ortho Clinical Diagnostics is unable to provide a definitive statement on the performance of its products with these devices. Confirm that your collection devices are compatible with this test.

Specimen Collection and Preparation

- Collect specimens using standard procedures.¹⁷
- Follow the instructions provided with your collection device for use and processing of the sample.¹⁸
- Samples should be thoroughly separated from all cellular material. Failure to do so may lead to an erroneous result.
- Serum and plasma (lithium heparin) are the recommended sample types. **Serum and lithium heparin samples should not be used interchangeably when collecting serial samples from the same patient.**
- The VITROS hs Troponin I test uses 60 µL of sample for each determination. This does not take account of the minimum fill volume of the chosen sample container. For details on minimum fill volume of sample cups or containers, refer to the operating instructions for your system.

Handling and Storage Conditions

- Handle samples in stoppered containers to avoid contamination and evaporation.
- Follow procedures within your laboratory to avoid cross contamination of patient specimens.
- The amount of time samples are on the system prior to analysis should be limited to avoid evaporation. Refer to the operating instructions for your system.
- Return to 2–8 °C (36–46 °F) as soon as possible after use, or load sufficient volume for a single determination.

- Serum and plasma samples may be stored for up to 8 hours at room temperature 15–30 °C (59–86 °F). Serum and plasma samples may be stored for up to 2 days at 2–8 °C (36–46 °F). Thoroughly mix samples by inversion and bring to 15–30 °C (59–86 °F) before use.
- Samples that will not be tested within the time frames outlined above should be stored frozen at -20 °C (-4 °F). Serum and plasma samples tested initially and up to 13 weeks storage at -20 °C (-4 °F) showed no differences in clinical performance. Serum and plasma samples may be subjected to up to two freeze-thaw cycles.
- After thawing, frozen samples must be mixed thoroughly by inversion and centrifuged before performing the assay. After centrifuging transfer the supernatant to a sample cup or secondary tube for testing to avoid any precipitate.

Testing Procedure

Materials Provided

- VITROS Immunodiagnostic Products hs Troponin I Reagent Pack
- VITROS Immunodiagnostic Products hs Troponin I Calibrators

Materials Required but Not Provided

- VITROS Immunodiagnostic Products Signal Reagent
- VITROS Immunodiagnostic Products Universal Wash Reagent
- VITROS Immunodiagnostic Products High Sample Diluent B
- VITROS Immunodiagnostic Products hs Troponin I Low Control
- VITROS Immunodiagnostic Products Reagent Pack Storage Box with desiccant
- Other quality control materials

Operating Instructions

Check the inventory regularly to aid the management of reagents and ensure that sufficient VITROS Signal Reagent, VITROS Universal Wash Reagent and calibrated reagent lots are available for the work planned. When performing panels of tests on a single sample, ensure that the sample volume is sufficient for the tests ordered. Ensure sufficient VITROS High Sample Diluent B Reagent Pack is loaded onto the system before processing samples. Refer to the VITROS High Sample Diluent B Reagent Pack instructions for use.

For detailed information refer to the operating instructions for your system.

Note: Do not use visibly damaged product.

Sample Dilution

Samples with concentrations greater than the measuring range may be automatically diluted on the system up to 5-fold (1 part sample with 4 parts diluent) by the VITROS Immunodiagnostic and VITROS Integrated Systems with the VITROS Immunodiagnostic Products High Sample Diluent B Reagent Pack prior to testing. Refer to the VITROS High Sample Diluent B Reagent Pack Instructions for Use.

Default Test Name

The default test name which will appear on patient reports is hs Troponin I. The default short name that will appear on the test selection menus and laboratory reports is hsTnI. These defaults may be reconfigured, if required. For detailed information refer to the operating instructions for your system.

Calibration

Calibration Procedure

- Calibration is lot specific; reagent packs and calibrators are linked by lot number. Reagent packs from the same lot may use the same calibration.
- A Master Calibration (a dose response curve covering the full calibration range) is established for each new reagent lot. Concentrations for the linked lot of calibrators are determined from the Master Calibration.
- Ensure that the Master Calibration for each new reagent lot is available on your system.
- Process calibrators in the same manner as samples. Calibration need not be programmed if bar code labels are used; load the calibrators in any order, calibration will be initiated automatically.
- When the calibrators are processed the signal expected for each calibrator is compared against the actual signal obtained. The Master Calibration is then rescaled to reflect the differences between the actual and expected signals. The validity of this calibration curve is assessed against a range of quality parameters, and if acceptable, it is stored for use with any reagent pack of that lot.

- The quality of calibration cannot be completely described by a single parameter. The calibration report should be used in conjunction with acceptable control values to determine the validity of the calibration.
- Recalibration is required after a pre-determined calibration interval, or when a different reagent lot is loaded.
- Calibration results are assessed against a range of quality parameters. Failure to meet any of the defined quality parameter ranges will be coded in the calibration report. For actions to be taken following a failed calibration refer to the operating instructions for your system.

Refer to the operating instructions for your system for detailed instructions on the calibration process.

When to Calibrate

- Calibrate when the reagent pack and calibrator lot changes.
- Calibrate every 28 days.
- After specified service procedures have been performed.
- If quality control results are consistently outside of your acceptable range.

For additional information on when to calibrate, refer to the operating instructions for your system.

Traceability of Calibration

Calibration of the VITROS hs Troponin I test is traceable to in-house reference calibrators, which have been value assigned to correlate to another commercially available test.

Calibration Model

A modified four-parameter logistic curve fit function is used to construct the Master Calibration. The calibration process rescales the Master Calibration to establish a valid stored curve for the VITROS Immunodiagnostic System and VITROS Integrated Systems.

Measuring (Reportable) Range

VITROS System	Measuring (Reportable) Range
ECi/ECiQ, 3600, 5600, XT 7600	1.50–30,000 ng/L (1.50–30,000 pg/mL)

Quality Control

Quality Control Material Selection

VITROS hs Troponin I Low Control and other controls containing suitable levels of cardiac troponin I are recommended for use with the VITROS Immunodiagnostic and VITROS Integrated Systems. The VITROS hs Troponin I Low Control contains 1 level troponin I (targeted near 99th percentile). The performance of commercial control fluids should be evaluated for compatibility with this test before they are used for quality control.

Control materials may show a difference when compared with other cardiac troponin I methods if they contain high concentrations of preservatives, stabilizers, or other nonphysiological additives, or otherwise depart from a true human sample matrix.

Appropriate quality control value ranges must be established for all quality control materials used with the VITROS hs Troponin I test.

Quality Control Procedure Recommendations

- Good laboratory practice requires that controls be processed to verify the performance of the test.
- Choose control levels that check the clinically relevant concentrations.
- To verify system performance, analyze control materials:
 - After calibration
 - If the system is turned off for more than 2 hours
 - After reloading reagent packs that have been removed from the MicroWell Supply and stored for later use
 - According to local regulations or at least once each day that the test is being performed
 - After specified service procedures are performed

If quality control procedures within your laboratory require more frequent use of controls, follow those procedures.

- Analyze quality control materials in the same manner as patient specimens.
- If control results fall outside your acceptable range, investigate the cause before deciding whether to report patient results.
- Refer to published guidelines for general quality control recommendations. ¹⁹

For more detailed information, refer to the operating instructions for your system.

Quality Control Material Preparation and Storage

Refer to the manufacturer's product literature for preparation, storage, and stability information.

Results

Results are automatically calculated by the VITROS Immunodiagnostic and VITROS Integrated Systems.

Reporting Units and Unit Conversion

Analyte results are quoted in units of ng/L or pg/mL. To configure the units, refer to the operating instructions for your system.

Conventional	Alternate
ng/L (pg/mL × 1)	pg/mL (ng/L × 1)

Limitations of the Procedure

Known Interferences

The VITROS hs Troponin I test was evaluated for interference consistent with CLSI document EP07²⁰. Commonly encountered substances were tested on three lots of reagents. The following compounds, when tested, caused the bias shown at the concentrations indicated.

Refer to "Specificity" for a list of compounds tested that did not show interference.

Interferent	Interferent Concentration		Analyte Conc [*]	% Bias ^{**}
			ng/L	
Acetylcysteine	15.0 mg/dL	920 µmol/L	14.83	17.2
			348.4	11.7
	6.83 mg/dL	419 µmol/L	13.98	10.0
	9.95 mg/dL	610 µmol/L	326.8	10.0
Bilirubin, unconjugated	40 mg/dL	684 µmol/L	9.22	-25.8
	10.4 mg/dL	178 µmol/L	9.53	-10.0
Cholesterol	400 mg/dL	10.3 mmol/L	8.40	-31.8
	315 mg/dL	8.15 mmol/L	9.43	-10.0
Dextran	2400 mg/dL	600 µmol/L	18.27	89.5
			320.4	34.7
	433 mg/dL	108 µmol/L	13.39	10.0
Fibrinogen	1000 mg/dL	29.4 µmol/L	8.32	-26.6
	490 mg/dL	14.4 µmol/L	9.13	-10.0
Hemoglobin	1000 mg/dL	10.0 g/L	213.4	-15.3
	608 mg/dL	6.08 g/L	259.6	-10.0
Sodium Azide	100 mg/dL	15,382 µmol/L	9.78	-19.1
	33.9 mg/dL	5215 µmol/L	11.43	-10.0
Streptokinase	150,000 U/dL		8.56	-34.0
			18.56	-94.1
	11,111 U/dL	11.38	-10.0	
Total Protein	4167 U/dL		226.4	-10.0
	15.0 g/dL		9.74	-28.1
	11.2 g/dL		274.2	-16.2
	14.5 g/dL		8.93	-10.0
14.5 g/dL		249.9	-10.0	

^{*} Average test concentration of replicate determinations using 3 different lots of reagent.

^{**} Estimate of the maximum difference observed as a percentage.

Note: These results are representative. The degree of interference at concentrations other than those listed might not be predictable from these results. Other interfering substances may be encountered in the patient population.

Other Limitations

- The VITROS hs Troponin I test has no high dose hook effect up to a concentration of 100,000 ng/L (100,000 pg/mL).
- Heterophilic antibodies in serum or plasma samples may cause interference in immunoassays.²¹ These antibodies may be present in blood samples from individuals regularly exposed to animals or who have been treated with animal serum products. Results that are inconsistent with clinical observations indicate the need for additional testing.

Expected Values and Interpretation of Results

The ESC/ACC/AHA/WHF Expert Consensus Document established that acute myocardial injury with clinical, ECG and/or imaging evidence of acute myocardial ischemia and rise and/or fall of cardiac troponin (cTn) with at least one value above the 99th percentile Upper Reference Limits (URLs) are criteria for defining acute myocardial infarction (AMI).⁵

Expected Values:

It is recommended that each laboratory establish its own upper reference limits and confirm the validity of the diagnostic cutoff for the population it serves.

The VITROS hs Troponin I 99th percentile URLs were established from matched serum and lithium heparin plasma of nine hundred fifty-two (952) apparently healthy adults, including 486 female and 466 male subjects. The results included variation from reagent lots, systems, operators, and reagent age. The subjects ranged in age from 22 to 91 years old, with fifty-nine percent of the subjects ≥ 50 years of age.

Subjects were excluded if they met any of the following exclusion criteria:

- History of kidney disease, diabetes, heart disease, cancer, lung disease, thyroid disease, or stroke
- High blood pressure, cholesterol or triglycerides
- Muscle or skeletal injury or surgery in the last 3 months
- Current smoker
- Pregnancy
- Additional exclusion criteria:
 - Hemoglobin A1c ≥ 6.5%
 - NT-proBNP > 125 pg/mL for subjects < 75 years of age or > 450 pg/mL for subjects ≥ 75 years of age
 - eGFR < 60 mL/min

The overall observed 99th percentile URL from the 952 lithium heparin plasma samples is 11 ng/L (90% CI of 8.4–14.8 ng/L). The overall observed 99th percentile from the 952 serum samples is 11 ng/L (90% CI of 7.8–16.3 ng/L). The 99th percentile URL values and respective 90% Confidence Intervals (CI), determined for lithium heparin plasma (females, males, and overall) and serum (females, males, and overall) using the non-parametric statistical method, are shown in the table below.

Sample Type	Gender	Number of Samples	99 th Percentile URL ng/L (pg/mL)	90% CI* ng/L (pg/mL)
Lithium Heparin Plasma	Female	486	9	4.8–19.0
	Male	466	13	9.7–22.9
	Overall	952	11	8.4–14.8
Serum	Female	486	9	4.4–21.4
	Male	466	12	8.8–23.7
	Overall	952	11	7.8–16.3

* CI = Confidence Interval

Alternate units: 1.00 ng/L = 1.00 pg/mL

For a troponin test to be labelled as high sensitivity it has to fulfil the two IFCC TF-CB established criteria²²:

- The test should have analytical imprecision (% CV) at the 99th percentile ≤ 10% and
- The test should measure troponin concentrations above the Limit of Detection (LoD) in ≥ 50% of healthy subjects

The VITROS hs Troponin I test met these criteria. The 10% CV Limit of Quantitation (LoQ) was measured to be 1.99 ng/L (1.99 pg/mL). Greater than 50% of the healthy subjects used to establish the 99th percentile had cTnI results above the claimed LoDs on the VITROS Immunodiagnostic and VITROS Integrated Systems.

Interpretation of Results:

Myocardial injury is defined by elevated cTn values above the 99th percentile and can occur in various cardiac and/or systemic conditions not associated with AMI. Conditions such as heart failure, myocarditis, structural heart disease, sepsis, kidney disease, pulmonary embolism, pulmonary hypertension, stroke, chemotherapeutic agents, critically ill patients, and strenuous exercise can result in acute or chronic myocardial injury. While in acute myocardial injury cTn levels are expected to rise and/or fall, in chronic myocardial injury such as structural heart disease and chronic kidney disease, the cTn levels are stable.

Due to the multitude of conditions associated with troponin elevations, serial measurements of cTn are recommended and detection of a rise and/or fall of cardiac troponin (cTn) values in the setting of acute myocardial ischemia is considered essential for AMI diagnosis.⁵

Results of this test should be used in conjunction with clinical presentation, other diagnostic tests, and in accordance with the appropriate clinical guidelines.

Performance Characteristics

Clinical Sensitivity and Specificity

A subset of serial serum samples obtained from the Advantageous Predictors of Acute Coronary Syndromes Evaluation (APACE) study, a well-characterized, ongoing prospective, international, multicenter sample collection in 12 sites in five European countries, was used in the clinical validation. See Trial Registration No: NCT00470587 at clinicaltrials.gov for more details. The enrolled subjects presented to the emergency department with symptoms of chest pain and angina pectoris. Peak symptoms had to have occurred within the last 12 hours. Subjects with kidney failure requiring dialysis were excluded. Information from subjects in the study was adjudicated by at least two independent cardiologists and a consensus diagnosis was obtained. Diagnosis of MI was conducted using the diagnostic criteria described in the ACC/ECS/AHA guidelines.^{23,24} Once the subjects were discharged they were contacted by telephone at 30 days, 6 months, and 12 months to determine major adverse cardiac events (cardiac-related mortality and myocardial infarction) and all-cause mortality.

The 1231 subjects in the study ranged in age from 18 years to 93 years (females 18-92 years, males 19-93 years). The population consisted of 420/1231 (34.12%) females and 811/1231 (65.88%) males. A final adjudicated diagnosis of MI was assigned to 158/1231 (12.84%) subjects and a final adjudicated diagnosis of no MI was assigned to 1073/1231 (87.16%) subjects. Of the 158 MI subjects, 41 (25.95%) were female and 117 (74.05%) were male, and of the 1073 no MI subjects, 379 (35.32%) were female and 694 (64.68%) were male. The prevalence of MI in females was 9.76% (41/420) and in males was 14.43% (117/811).

Serial serum samples were tested and analyzed for each time point within and across gender. Receiver Operating Characteristic (ROC) curves²⁵ were generated. The diagnostic accuracy measurements of the VITROS hs Troponin I test, as quantified by areas under the ROC curves (AUC),²⁶ for each time point within and across gender ranged from 0.94 to 0.97.

The clinical performance [clinical sensitivity, clinical specificity, negative predictive value (NPV), positive predictive value (PPV), negative likelihood ratio (LR-) and positive likelihood ratio (LR+)] and the two-tailed 95% exact confidence intervals (CI) of the VITROS hs Troponin I test versus adjudicated diagnosis was determined using the overall and gender-specific 99th percentile cutoffs and are summarized in the following tables:

All Subjects using 11 ng/L Cutoff (Overall 99th Percentile - Serum)

Time Point	N	Sensitivity (%) (n/N)	95% CI (%)	Specificity (%) (n/N)	95% CI (%)	NPV (%) (n/N)	95% CI (%)	PPV (%) (n/N)	95% CI (%)
Baseline	1231	81.65 (129/158)	74.72–87.35	91.99 (987/1073)	90.20–93.54	97.15 (987/1016)	95.93–98.08	60.00 (129/215)	53.12–66.60
1 hour	1039	93.75 (120/128)	88.06–97.26	90.67 (826/911)	88.59–92.48	99.04 (826/834)	98.12–99.58	58.54 (120/205)	51.47–65.36
2 hours	869	95.79 (91/95)	89.57–98.84	89.53 (693/774)	87.16–91.60	99.43 (693/697)	98.54–99.84	52.91 (91/172)	45.16–60.55
3 hours	442	94.44 (51/54)	84.61–98.84	87.63 (340/388)	83.93–90.74	99.13 (340/343)	97.47–99.82	51.52 (51/99)	41.25–61.68

Time Point	N	LR-	95% CI	LR+	95% CI
Baseline	1231	0.20	0.14–0.28	10.19	8.21–12.64
1 hour	1039	0.07	0.04–0.13	10.05	8.17–12.36
2 hours	869	0.05	0.02–0.12	9.15	7.42–11.30
3 hours	442	0.06	0.02–0.19	7.63	5.81–10.03

All Subjects using 9 ng/L Cutoff for Females and 12 ng/L for Males (Gender-Specific 99th Percentiles - Serum)

Time Point	N	Sensitivity (%) (n/N)	95% CI (%)	Specificity (%) (n/N)	95% CI (%)	NPV (%) (n/N)	95% CI (%)	PPV (%) (n/N)	95% CI (%)
Baseline	1231	82.28 (130/158)	75.42–87.89	91.71 (984/1073)	89.89–93.29	97.23 (984/1012)	96.03–98.15	59.36 (130/219)	52.54–65.93
1 hour	1039	93.75 (120/128)	88.06–97.26	90.67 (826/911)	88.59–92.48	99.04 (826/834)	98.12–99.58	58.54 (120/205)	51.47–65.36
2 hours	869	93.68 (89/95)	86.76–97.65	89.66 (694/774)	87.30–91.72	99.14 (694/700)	98.14–99.68	52.66 (89/169)	44.85–60.38
3 hours	442	92.59 (50/54)	82.11–97.94	87.89 (341/388)	84.22–90.96	98.84 (341/345)	97.06–99.68	51.55 (50/97)	41.18–61.82

Time Point	N	LR-	95% CI	LR+	95% CI
Baseline	1231	0.19	0.14–0.27	9.92	8.03–12.26
1 hour	1039	0.07	0.04–0.13	10.05	8.17–12.36
2 hours	869	0.07	0.03–0.15	9.06	7.32–11.23
3 hours	442	0.08	0.03–0.22	7.64	5.79–10.10

Females using 9 ng/L Cutoff (Female 99th Percentile Cutoff-Serum)

Time Point	N	Sensitivity (%) (n/N)	95% CI (%)	Specificity (%) (n/N)	95% CI (%)	NPV (%) (n/N)	95% CI (%)	PPV (%) (n/N)	95% CI (%)
Baseline	420	85.37 (35/41)	70.83–94.43	92.88 (352/379)	89.80–95.25	98.32 (352/358)	96.39–99.38	56.45 (35/62)	43.26–69.01
1 hour	346	94.12 (32/34)	80.32–99.28	91.99 (287/312)	88.40–94.75	99.31 (287/289)	97.52–99.92	56.14 (32/57)	42.36–69.26
2 hours	291	91.67 (22/24)	73.00–98.97	91.01 (243/267)	86.92–94.16	99.18 (243/245)	97.08–99.90	47.83 (22/46)	32.89–63.05
3 hours	133	92.86 (13/14)	66.13–99.82	88.24 (105/119)	81.05–93.42	99.06 (105/106)	94.86–99.98	48.15 (13/27)	28.67–68.05

Time Point	N	LR-	95% CI	LR+	95% CI
Baseline	420	0.16	0.08–0.33	11.98	8.15–17.61
1 hour	346	0.06	0.02–0.25	11.75	7.99–17.27
2 hours	291	0.09	0.02–0.35	10.20	6.83–15.22
3 hours	133	0.08	0.01–0.54	7.89	4.73–13.18

Males using 12 ng/L Cutoff (Male 99th Percentile Cutoff-Serum)

Time Point	N	Sensitivity (%) (n/N)	95% CI (%)	Specificity (%) (n/N)	95% CI (%)	NPV (%) (n/N)	95% CI (%)	PPV (%) (n/N)	95% CI (%)
Baseline	811	81.20 (95/117)	72.93–87.82	91.07 (632/694)	88.69–93.08	96.64 (632/654)	94.95–97.88	60.51 (95/157)	52.41–68.21
1 hour	693	93.62 (88/94)	86.62–97.62	89.98 (539/599)	87.29–92.27	98.90 (539/545)	97.62–99.59	59.46 (88/148)	51.09–67.44
2 hours	578	94.37 (67/71)	86.20–98.44	88.95 (451/507)	85.90–91.55	99.12 (451/455)	97.76–99.76	54.47 (67/123)	45.25–63.47
3 hours	309	92.50 (37/40)	79.61–98.43	87.73 (236/269)	83.20–91.40	98.74 (236/239)	96.38–99.74	52.86 (37/70)	40.55–64.91

Time Point	N	LR-	95% CI	LR+	95% CI
Baseline	811	0.21	0.14–0.30	9.09	7.06–11.71
1 hour	693	0.07	0.03–0.15	9.35	7.31–11.95
2 hours	578	0.06	0.02–0.16	8.54	6.63–11.01
3 hours	309	0.09	0.03–0.25	7.54	5.41–10.50

The cutoffs and clinical performance information above should only be used as a guide. It is recommended that each laboratory establish its own upper reference limits and confirm the validity of the diagnostic cutoff for the population it serves.

Results of this test should be used in conjunction with clinical presentation, other diagnostic tests, and in accordance with the appropriate clinical guidelines.

Risk Stratification

The VITROS hs Troponin I test was evaluated for its prognostic value as an aid in the assessment of cardiac-related mortality/MI and all-cause mortality over 30 days (d), 6 months (m) and 12 months. Of the 1231 subjects enrolled in the study, 1226 were followed up. The cumulative event free rates (CEFR) (Kaplan-Meier analysis) over 30 days, 6 months and 12 months for subjects classified according to the overall and gender-specific 99th percentile cutoffs are summarized in the following tables:

Prognosis for Cardiac-Related Mortality and MI

VITROS hs Troponin I Test	Follow up period	≤ Cutoff				> Cutoff				Log Rank p-value
		Subjects at Risk (N)	Events (N)	Censored (N)	CEFR (%)	Subjects at Risk (N)	Events (N)	Censored (N)	CEFR (%)	
Overall 99 th Percentile 11 ng/L	30 d	966	3	963	99.69	260	5	255	98.08	0.0041
	6 m	966	9	957	99.07	260	11	249	95.77	0.0002
	12 m	966	14	952	98.55	260	14	246	94.62	0.0002
Gender-Specific 99 th Percentile 9 ng/L female 12 ng/L male	30 d	968	4	964	99.59	258	4	254	98.45	0.0439
	6 m	968	10	958	98.97	258	10	248	96.12	0.0014
	12 m	968	15	953	98.45	258	13	245	94.96	0.0010

Prognosis for All-Cause Mortality

VITROS hs Troponin I Test	Follow up period	≤ Cutoff				> Cutoff				Log Rank p-value
		Subjects at Risk (N)	Deaths (N)	Censored (N)	CEFR (%)	Subjects at Risk (N)	Deaths (N)	Censored (N)	CEFR (%)	
Overall 99 th Percentile 11 ng/L	30 d	966	1	965	99.90	260	4	256	98.46	0.0013
	6 m	966	9	957	99.07	260	8	252	96.92	0.0083
	12 m	966	15	951	98.45	260	12	248	95.38	0.0030
Gender-Specific 99 th Percentile 9 ng/L female 12 ng/L male	30 d	968	1	967	99.90	258	4	254	98.45	0.0012
	6 m	968	9	959	99.07	258	8	250	96.90	0.0079
	12 m	968	14	954	98.55	258	13	245	94.96	0.0006

The cohort of subjects with at least one VITROS hs Troponin I test value above the overall and gender-specific cutoffs had significantly lower cumulative event free rates for cardiac-related mortality/MI and all-cause mortality at 30 days, 6 months and 12 months than the cohort of subjects with all VITROS hs Troponin I test values less than or equal to the overall and gender-specific cutoffs.

Results of this test should be used in conjunction with clinical presentation, other diagnostic tests, and in accordance with the appropriate clinical guidelines.

Limit of Detection

The Limit of Detection (LoD) for the VITROS hs Troponin I test was designed to be < 3.00 ng/L in accordance with CLSI document EP17²⁷. The LoD for the VITROS hs Troponin I test is between 0.39 and 0.86 ng/L, based on 500 determinations with five low-level samples. The Limit of Blank (LoB) is between 0.14 and 0.51 ng/L, based on 400 determinations of four blank samples. Representative performance data are shown below.

Limit of Blank and Limit of Detection

System	LoB		LoD	
	ng/L	pg/mL	ng/L	pg/mL
ECi/ECiQ	0.51	0.51	0.86	0.86
3600	0.19	0.19	0.39	0.39
5600	0.23	0.23	0.43	0.43
XT 7600	0.14	0.14	0.42	0.42

Limit of Quantitation

The Limit of Quantitation (LoQ) for the VITROS hs Troponin I test is 1.23 ng/L, as determined based on CLSI document EP17²⁷, from 75 determinations with five LoQ pools and a precision goal of 20% CV. The LoQ was determined to be the lowest concentration at which precision and accuracy design requirements are still met and within the linear range of the test for the VITROS ECi/ECiQ and 3600 Immunodiagnostic Systems and the VITROS 5600 and XT 7600 Integrated Systems.

Limit of Quantitation

LoQ (20% CV)	
ng/L	pg/mL
1.23	1.23

Accuracy

Accuracy was evaluated consistent with CLSI document EP09.²⁸ The table shows the results of endogenous serum samples analyzed on the VITROS ECi/ECiQ Immunodiagnostic System compared with those analyzed using the VITROS 3600 Immunodiagnostic System, the VITROS 5600 Integrated System, and the VITROS XT 7600 Integrated System. The relationship between the VITROS systems was determined by Weighted Deming regression and Pearson correlation.

VITROS System	n	Slope	Correlation Coefficient	Conventional Units (ng/L)*	
				Range of Samples	Intercept
3600 vs. ECi/ECiQ	112	1.01	1.00	1.20–26340	-0.09
5600 vs. ECi/ECiQ	112	0.99	1.00	1.20–26340	-0.17
XT 7600 vs. ECi/ECiQ	112	1.03	1.00	1.11–26890	-0.24
5600 vs. 3600	112	0.98	1.00	1.20–26170	-0.08
XT 7600 vs. 3600	112	1.02	1.00	1.11–26890	-0.15
XT 7600 vs. 5600	112	1.04	1.00	1.11–26890	-0.07

* Alternate units: 1.00 ng/L = 1.00 pg/mL

Precision

Precision was evaluated consistent with CLSI document EP05.²⁹ Two replicates each of eight patient pools (four human serum and four human plasma) and two controls were tested on two separate occasions per day on at least 20 different test days. The experiment was performed using three reagent lots on one VITROS ECi/ECiQ Immunodiagnostic System, one VITROS 3600 Immunodiagnostic System, one VITROS 5600 Integrated System, and one VITROS XT 7600 Integrated System. The test was designed to have within-laboratory precision of $\leq 20\%$ at or near the LoQ, $\leq 10\%$ at the 99th percentile, and $< 7\%$ between the 99th percentile and 30,000 ng/L. Representative performance data are shown below.

VITROS System	Conc. Units - ng/L							No. Observ.	No. Days
	Mean hs Troponin I Conc.	Within-run*		Within-calibration**		Within-lab***			
		SD	CV (%)	SD	CV (%)	SD	CV (%)		
ECi/ECiQ	2.97	0.217	7.1	0.286	9.4	0.403	13.9	80	20
	10.75	0.359	3.3	0.612	5.6	0.785	7.4	80	20
	63.22	1.299	2.1	2.768	4.4	3.074	4.9	80	20
	303.7	5.73	1.9	8.57	2.8	8.82	2.9	80	20
	12000	316.1	2.7	454.7	3.9	540.7	4.4	80	20
	16380	461.3	2.9	794.7	5.0	854.5	5.0	80	20
	2.19	0.202	9.0	0.284	12.7	0.382	17.9	80	20
	6.68	0.276	4.1	0.487	7.2	0.612	9.3	80	20
	56.86	1.759	3.1	2.739	4.8	2.951	5.2	80	20
	298.0	6.71	2.3	10.62	3.6	11.25	3.7	80	20

VITROS System	Conc. Units - ng/L							No. Observ.	No. Days
	Mean hs Troponin I Conc.	Within-run*		Within-calibration**		Within-lab***			
		SD	CV (%)	SD	CV (%)	SD	CV (%)		
3600	2.53	0.055	2.3	0.167	6.9	0.258	9.8	80	20
	10.24	0.239	2.4	0.528	5.3	0.702	6.7	80	20
	67.55	0.779	1.2	2.558	3.8	3.097	4.5	80	20
	325.5	4.28	1.3	10.02	3.1	12.79	3.9	80	20
	14220	407.8	2.9	548.7	3.9	659.6	4.6	80	20
	18860	343.5	1.8	794.4	4.2	992.0	5.3	80	20
	1.84	0.064	3.7	0.116	6.6	0.192	9.9	80	20
	5.59	0.092	1.7	0.317	5.8	0.436	7.6	80	20
	65.60	0.810	1.3	2.518	3.9	3.300	5.0	80	20
5600	323.2	6.15	1.9	12.9	4.0	15.56	4.8	80	20
	2.44	0.076	3.3	0.152	6.6	0.229	8.9	80	20
	10.13	0.191	2.0	0.395	4.1	0.640	6.1	80	20
	66.46	1.142	1.8	2.625	4.0	3.051	4.5	80	20
	308.5	4.72	1.6	9.00	3.0	9.02	2.9	80	20
	14090	262.7	1.9	415.8	3.0	764.8	5.4	80	20
	18460	236.3	1.3	565.0	3.1	544.8	2.9	80	20
	1.71	0.119	7.4	0.151	9.4	0.206	11.3	80	20
	5.48	0.141	2.7	0.254	4.9	0.422	7.4	80	20
XT 7600	65.27	1.143	1.8	1.907	3.0	2.602	3.9	80	20
	309.2	3.42	1.1	9.15	3.0	9.56	3.0	80	20
	2.55	0.070	2.9	0.125	5.2	0.248	9.3	80	20
	10.30	0.271	2.8	0.464	4.7	0.920	8.5	80	20
	66.92	1.267	2.0	2.567	4.0	4.306	6.2	80	20
	310.2	4.02	1.3	10.60	3.5	14.59	4.6	80	20
	13750	260.8	1.9	552.7	4.1	616.8	4.4	80	20
	18360	461.1	2.5	776.2	4.3	844.6	4.5	80	20
	1.82	0.090	5.2	0.119	6.9	0.174	9.1	80	20
5.57	0.114	2.2	0.267	5.1	0.531	9.1	80	20	
65.50	0.914	1.4	2.407	3.8	4.064	6.0	80	20	
307.7	3.97	1.3	10.60	3.5	14.71	4.7	80	20	

* Within-run (repeatability). Between Duplicate precision averaged over all runs

** Within-calibration. Total precision with weighted components of within-run, between-run, and between-day variation.

*** Within-lab. A measure of the effect of recalibration on total precision, calculated within reagent lot, using data from at least 5 calibrations

Specificity

Substances that Do Not Interfere

The VITROS hs Troponin I test was evaluated for interference consistent with CLSI document EP07²⁰. Of the compounds tested, none was found to cause a bias of >10% with the test at the concentrations indicated at nominal troponin I concentrations of 10.00 ng/L and 300.0 ng/L.

Compound	Concentration	
Acetaminophen	156 µg/mL	1030 µmol/L
Adrenaline (epinephrine)	20 µg/dL	1.09 µmol/L
Allopurinol	6.0 mg/dL	441 µmol/L
Alprazolam	0.0258 mg/dL	0.835 µmol/L
Ambroxol	63 µg/dL	1.52 µmol/L
Amlodipine besylate	0.0104 mg/dL	0.183 µmol/L
Amoxicillin	5.40 mg/dL	148 µmol/L
Ascorbic acid	5.25 mg/dL	298 µmol/L
Atorvastatin calcium	0.162 mg/dL	1.34 µmol/L
Benazepril HCl	0.044 mg/dL	0.954 µmol/L

Compound	Concentration	
Bilirubin, conjugated	40 mg/dL	475 µmol/L
Biotin	3510 ng/mL	14,300 nmol/L
Bivalirudin	2.18 mg/dL	10.0 µmol/L
Caffeine	10.8 mg/dL	556 µmol/L
Carvedilol	43.2 µg/dL	1.06 µmol/L
Cefoxitin sodium	697 mg/dL	15,500 µmol/L
Ceftriaxone disodium hemi (heptahydrate)	100 mg/dL	1510 µmol/L
Cephalexin sodium	13.4 mg/dL	363 µmol/L
Cinnarizine	108 µg/dL	2.93 µmol/L
Clopidogrel	2.4 µg/dL	0.057 µmol/L
Cocaine	0.6 mg/dL	17.7 µmol/L
Cotinine	0.24 mg/dL	4.89 µmol/L
Cyclosporin A	0.18 mg/dL	1.50 µmol/L
Digoxin	0.0039 mg/dL	0.0499 µmol/L
Diphenhydramine	0.0774 mg/dL	3.03 µmol/L
Dopamine	0.0621 mg/dL	4.06 µmol/L
Enalaprilat	0.0819 mg/dL	2.35 µmol/L
Enoxaparin (LMWH)	360 U/dL	N/A
Eptifibatide	1.44 mg/dL	17.3 µmol/L
Erythromycin	13.8 mg/dL	188 µmol/L
Ethanol	600 mg/dL	130,000 µmol/L
Fondaparinux	0.39 mg/dL	2.26 µmol/L
Furosemide	1.59 mg/dL	48.1 µmol/L
HAMA (human anti-mouse antibodies)	800 µg/L	0.005 µmol/L
Heparin (Sodium), UFH	330 U/dL	N/A
Ibuprofen	21.9 mg/dL	1060 µmol/L
Insulin	3.12 µg/dL	0.005 µmol/L
L-dopa (Levodopa)	0.75 mg/dL	38.0 µmol/L
Levothyroxine	0.0429 mg/dL	0.552 µmol/L
Lidocaine	1.5 mg/dL	64.0 µmol/L
Methylprednisolone	0.783 mg/dL	20.9 µmol/L
Metronidazole	12.3 mg/dL	719 µmol/L
Naproxen sodium	39.3 mg/dL	1560 µmol/L
Nifedipine	0.0588 mg/dL	1.70 µmol/L
Nitrofurantoin	0.213 mg/dL	8.94 µmol/L
Nitroglycerin (Nitrostat)	1.2 µg/dL	0.053 µmol/L
Omeprazole	0.84 mg/dL	24.3 µmol/L
Oxycodone	0.0324 mg/dL	1.03 µmol/L
Oxytetracycline	1.2 mg/dL	24.2 µmol/L
Phenytoin	6.0 mg/dL	238 µmol/L
Propranolol HCl	0.115 mg/dL	3.88 µmol/L
Pseudoephedrine	0.330 mg/dL	20.0 µmol/L
Quinidine	1.5 mg/dL	46.2 µmol/L
Rheumatoid Factor	900 IU/mL	N/A
Rifampicin (Rifampin)	4.8 mg/dL	58.3 µmol/L
Rivaroxaban	0.270 mg/dL	6.19 µmol/L
Salicylic acid	2.86 mg/dL	207 µmol/L
Spironolactone	0.0555 mg/dL	1.33 µmol/L
Theophylline	6.0 mg/dL	333 µmol/L
TPA (Alteplase)	1.2 mg/dL	0.171 µmol/L

Compound	Concentration	
	mg/dL	μmol/L
Triglyceride	1500 mg/dL	16.9 mmol/L
Vancomycin hydrochloride	12.3 mg/dL	82.8 μmol/L
Verapamil	0.16 mg/dL	3.51 μmol/L
Vorapaxar	36 μg/dL	0.731 μmol/L
Warfarin sodium	8.0 mg/dL	259 μmol/L

Cross-Reactivity

The cross-reactivity of the VITROS hs Troponin I test was evaluated by adding the following substances to one human serum sample pool containing no Troponin I.

Cross-Reactant	Cross Reactant Concentration		Mean Result of Control Pool	Mean Result of Cross-Reactant Pool	% Cross-Reactivity
	ng/mL	ng/L	ng/L	ng/L	
Actin (from Rabbit Muscle)	100	100,000	*	*	*
Cardiac Troponin C (Recombinant)	100	100,000	*	*	*
Cardiac Troponin T (Recombinant)	100	100,000	*	*	*
CK-MB (Recombinant)	100	100,000	*	*	*
Myoglobin (Recombinant)	1000	1000,000	*	*	*
Myosin (Recombinant)	100	100,000	*	*	*
Skeletal Troponin I	100	100,000	*	*	*
Tropomyosin (from porcine muscle)	200	200,000	*	*	*

* Not Detectable. Concentration was below the measuring range of the test, 1.50–30,000 ng/L.

The cross-reactivity of the VITROS hs Troponin I test was evaluated by adding the following substances to one human serum sample pool containing Troponin I at a concentration of 10.00 ng/L.

Cross-Reactant	Cross Reactant Concentration		Mean Result of Control Pool	Mean Result of Cross-Reactant Pool	% Cross-Reactivity
	ng/mL	ng/L	ng/L	ng/L	
Actin (from Rabbit Muscle)	100	100,000	16.98	12.23	0.0
Cardiac Troponin C (Recombinant)	100	100,000	16.98	13.57	0.0
Cardiac Troponin T (Recombinant)	100	100,000	16.98	12.16	0.0
CK-MB (Recombinant)	100	100,000	16.98	16.94	0.0
Myoglobin (Recombinant)	1000	1000,000	16.98	17.07	0.0
Myosin (Recombinant)	100	100,000	16.98	12.59	0.0
Skeletal Troponin I	100	100,000	16.98	12.66	0.0
Tropomyosin (from porcine muscle)	200	200,000	16.98	16.78	0.0

Cross-reactivity was expressed as the mean result obtained for the cross-reactant pool minus the mean result obtained for the control sample divided by the cross-reactant concentration in percentage term.




































$$\% \text{ Cross-reactivity} = \frac{(\text{Mean Troponin I Result Cross-reactant Pool}) - (\text{Mean Troponin I Result Control Sample})}{\text{Concentration of Cross-reactant}} \times 100$$

References

1. Tobacman LS. Thin Filament-Mediated Regulation of Cardiac Contraction. *Annu Rev Physiol*. 58:447-481; 1996.
2. Wilkinson JM, Grand JA. Comparison of Amino Acid Sequence of Tnl from Different Striated Muscles. *Nature*. 271:31-35;1978.
3. Bodor GS. Cardiac Tnl: A Highly Specific Biochemical Marker for Myocardial Infarction. *J Clin Immunoassay*. 17(1): 40-44; 1994.
4. Vallins WJ et al. Molecular Cloning of Human cTnl Using Polymerase Chain Reaction. *FEBS Letts*. 270(1,2):57-61; 1990.
5. Thygesen K, Alpert JS, Jaffe AS, et al. Fourth universal definition of myocardial infarction. *J Am Coll Cardiol* 72(18): 2231-2264; 2018.
6. Park K. C. et al., Cardiac troponins: from myocardial infarction to chronic disease *Cardiovascular Research* 113, 1708–1718; 2017.
7. Goodman SG, et al., GRACE Investigators. The diagnostic and prognostic impact of the redefinition of acute myocardial infarction: Lessons from the Global Registry of Acute Coronary Events (GRACE). *Am Heart J*. 151:654–60; 2006.
8. Heidenreich PA, Alloggiamento T, Melsop K, et al. The prognostic value of troponin in patients with non-ST elevation acute coronary syndromes: a meta analysis. *J Am Coll Cardiol*. 38:478-485; 2001.
9. Antman EM, Tanasijevic MJ, Thompson B, et al. Cardiac-Specific Troponin I Levels to Predict the Risk of Mortality in Patients with Acute Coronary Syndromes. *N Eng J Med*. 335:1342-1349; 1996.
10. Ottani F, Galvani M, Nicolini A, et al. Elevated cardiac troponin levels predict the risk of adverse outcome in patients with acute coronary syndrome. *Am Heart J*. 40:917-927; 2000.
11. Newby LK, Ohman EM, Biomarkers: Troponin testing--risk stratification to stratified medicine. *Nat Rev Cardiol*. Nov; 12(11):625-6; 2015.
12. Larue C et al. Cardiac-Specific Immunoenzymometric Assay of Tnl in the Early Phase of AMI. *Clin Chem*. 39(6): 972-979; 1993.
13. Roffi M et al., 2015 ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation. *European Heart Journal* 37, 267–315; 2016.
14. Summers M et al. Luminogenic Reagent Using 3-Chloro 4-Hydroxy Acetanilide to Enhance Peroxidase/Luminol Chemiluminescence. *Clin Chem*. 41:S73; 1995.
15. CLSI *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI document M29-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
16. Calam RR. Specimen Processing Separator Gels: An Update. *J Clin Immunoassay*. 11:86–90; 1988.
17. CLSI. *Collection of Diagnostic Venous Blood Specimens. 7th ed*. CLSI standard GP41. Wayne, PA: Clinical and Laboratory Standards Institute; 2017.
18. CLSI. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. CLSI document GP39-A6. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
19. CLSI. *Statistical Quality Control for Quantitative Measurement Procedures: Principles and Definitions; Approved Guideline - Fourth Edition*. CLSI guideline C24, Wayne, PA: Clinical and Laboratory Standards Institute; 2016.
20. CLSI. *Interference Testing in Clinical Chemistry. 3rd ed*. CLSI guideline EP07. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
21. Levinson SS. The Nature of Heterophilic Antibodies and Their Role in Immunoassay Interference. *J Clin Immunoassay*. 15:108-115; 1992.
22. Apple FS, Jaffe AS, Collinson P, et al. IFCC educational materials on selected analytical and clinical applications of high sensitivity cardiac troponin assays. *Clin Biochem* 48(4-5):201-3; 2015.
23. Thygesen K, Alpert JS, White HD, et al. Universal Definition of Myocardial Infarction. *Circulation*. 116:2634-2653; 2007.
24. Thygesen K, et al. Third Universal Definition of Myocardial Infarction. *Circulation*. 126:2020-2035; 2012.
25. Zweig MH, Campbell G. Receiver-Operating Characteristic (ROC) Plots: A Fundamental Evaluation Tool in Clinical Medicine. *Clin Chem*. 39(4):561-577; 1993.
26. Bamber D. The Area above the Ordinal Dominance Graph and the Area below the Receiver Operating Characteristic Graph. *Journal of Mathematical Psychology*. 12:387-415; 1975.
27. CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
28. CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples. 3rd ed*. CLSI guideline EP09. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
29. CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.

Glossary of Symbols

The following symbols may have been used in the labeling of this product.

	Do Not Reuse		Upper Limit of Temperature		Range
	Use by or Expiration Date (Year-Month-Day)		Lower Limit of Temperature		Range of Means
	Batch Code or Lot Number		Temperature Limitation		Midpoint
	Serial Number		Consult Instructions for Use		Revised
	Catalog Number or Product Code		Attention: The Instructions for Use (IFU) has been updated		Supersedes
	Caution		For use in Slide Supply 1		Contains Sufficient for "n" Tests
	Keep Dry (Protect from Moisture/Humidity)		For use in Slide Supply 2		<i>in vitro</i> Diagnostic Medical Device
	Manufacturer		SI Units		Der Grüne Punkt (the Green Dot). Manufacturer follows certain packaging material waste disposal management regulations
	Date of Manufacture		Conventional Units		Estimated within-lab SD
	Authorized Representative in the European Community		Value		Serious Health Hazards
	Corrosive		Flammable		Environmental or Aquatic Toxicity
	Health Hazards		Acute Toxicity		

Revision History

Date of Revision	Version	Description of Technical Changes*
2020-04-22	3.0	Calibrator Storage and Preparation: changed -20 °C to ≤-20 °C

*The change bars indicate the position of a technical amendment to the text with respect to the previous version of the document.

When this Instructions For Use is replaced, sign and date below and retain as specified by local regulations or laboratory policies, as appropriate.

Signature

Obsolete Date

Conditions of supply: all supplies are made subject to the standard terms and conditions of Ortho Clinical Diagnostics or its distributors. Copies of these are available on request.



Ortho-Clinical Diagnostics
1500 Boulevard Sébastien Brant
B.P. 30335
67411 Illkirch
CEDEX, France



Ortho-Clinical Diagnostics
Felindre Meadows
Pencoed
Bridgend
CF35 5PZ
United Kingdom

VITROS is a trademark of Ortho Clinical Diagnostics
The third party trademarks used herein are trademarks of their
respective owners.
© Ortho Clinical Diagnostics, 2019–2020