

INSTRUCTIONS FOR USE

Prog

VITROS Immunodiagnostic Products
Progesterone Reagent Pack

REF 193 8570

VITROS Immunodiagnostic Products
Progesterone Calibrators

REF 825 4898

Rx ONLY

Intended Use

For *in vitro* diagnostic use only.

VITROS Immunodiagnostic Products Progesterone Reagent Pack

For the quantitative measurement of progesterone in human serum and plasma (EDTA or heparin) using the VITROS ECI/ECiQ/3600 Immunodiagnostic Systems and the VITROS 5600/XT 7600 Integrated Systems.

VITROS Immunodiagnostic Products Progesterone Calibrators

For use in the calibration of the VITROS ECI/ECiQ/3600 Immunodiagnostic Systems and the VITROS 5600/XT 7600 Integrated Systems for the quantitative measurement of progesterone in human serum and plasma (EDTA or heparin).

Summary and Explanation of the Test

The main sites of progesterone production are the adrenal cortex, ovaries, and corpus luteum following ovulation and the placenta by the twelfth week of pregnancy.¹⁻³ Circulating progesterone is bound to several serum proteins including albumin and corticosteroid binding globulin. The physiologically active free hormone represents approximately 3% of the total progesterone concentration.⁴ Measurement of serum progesterone is useful in the investigation of ovarian function where disorders of ovulation are responsible for infertility in 15–20% of patients, and for predicting ovulation in induced cycles, where concentrations are generally higher than normal.⁵ Properly timed measurements of progesterone can be used in the diagnosis of patients with recurrent and threatened abortion in the first ten weeks of gestation.^{6,7} Corpus luteum dysfunction is indicated by lower than normal progesterone concentrations.

Principles of the Procedure

A competitive immunoassay technique is used. Progesterone present in the sample competes with a horseradish peroxidase (HRP)-labeled progesterone for a limited number of binding sites on a biotinylated rabbit anti-progesterone antibody presented in the liquid phase. The effects of binding proteins are eliminated by use of an appropriate blocking agent. The biotinylated rabbit anti-progesterone antibody is captured by streptavidin coated on the wells. Unbound materials are removed by washing.


The bound HRP conjugate is measured by a luminescent reaction.⁸ A reagent containing luminogenic substrates (a luminol derivative and a peracid salt) and an electron transfer agent, is added to the wells. The HRP in the bound conjugate catalyzes the oxidation of the luminol derivative, producing light. The electron transfer agent (a substituted acetanilide) increases the level of light produced and prolongs its emission. The light signals are read by the system. The amount of HRP conjugate bound is indirectly proportional to the concentration of progesterone present.

Test Type	System *	Incubation Time	Time to first result	Test Temperature	Reaction Sample Volume
Competitive immunoassay	ECI/ECiQ, 3600, 5600, XT 7600	16 minutes	24 minutes	37 °C	25 µL

* Not all products and systems are available in all countries.

Reaction Scheme

Streptavidin
Coated Well



Warnings and Precautions

WARNING: Potentially Infectious Material

Human blood products provided as components of the VITROS Progesterone Calibrators have been obtained from donors who were tested individually and who were found to be negative for hepatitis B surface antigen, and for antibodies to human immunodeficiency virus (HIV 1+2) and hepatitis C virus (HCV), using approved methods (enzyme immunoassays). Treat as if capable of transmitting infection.

Use caution when handling material of human origin. Consider all samples potentially infectious. No test method can offer complete assurance that hepatitis B virus, HCV, HIV 1+2 or other infectious agents are absent. Handle, use, store and dispose of solid and liquid waste from samples and test components, in accordance with procedures defined by appropriate national biohazard safety guideline or regulation (e.g. CLSI document M29).⁹

WARNING: Contains Kathon or ProClin 200 (CAS 55965-84-9)¹⁰

The VITROS Progesterone Reagent Pack and VITROS Progesterone Calibrators contain 1.0% Kathon or ProClin 200. H317: May cause an allergic skin reaction. P280: Wear protective gloves/protective clothing/eye protection/face protection. P302 + P352: IF ON SKIN: Wash with plenty of soap and water. P333 + P313: If skin irritation or rash occurs: Get medical advice/attention. P363: Wash contaminated clothing before reuse.

Refer to www.Orthoclinicaldiagnostics.com for the Safety Data Sheets and for Ortho contact information.

WARNING



Reagents

Reagent Pack Contents

- 1 reagent pack containing:
- 100 coated wells (streptavidin, bacterial; binds ≥ 3 ng biotin/well)
 - 10.0 mL conjugate reagent (HRP-progesterone, 70 ng/mL) in buffer with bovine serum albumin and antimicrobial agent
 - 10.0 mL biotinylated antibody reagent (biotin-rabbit polyclonal anti-progesterone, binds ≥ 70 pmol progesterone/mL) in buffer with bovine serum albumin, bovine gamma globulin and antimicrobial agent

Reagent Pack Handling

- The reagent pack is supplied ready for use.

- The reagent pack contains homogeneous liquid reagents that do not require shaking or mixing prior to loading onto the system.
- Handle the reagent pack with care. Avoid the following:
 - allowing condensation to form on the pack
 - causing reagents to foam
 - agitation of the pack

Reagent Pack Storage and Preparation

Reagent	Storage Condition		Stability
Unopened	Refrigerated	2–8 °C (36–46 °F)	expiration date
Opened	3600/5600/XT 7600: On system	System turned on	≤7 weeks
	ECi/ECiQ: On system	System turned on	≤12 weeks
Opened	3600/5600/XT 7600: Refrigerated	2–8 °C (36–46 °F)	≤7 weeks
	ECi/ECiQ: Refrigerated	2–8 °C (36–46 °F)	≤12 weeks

- The VITROS Progesterone Reagent Pack is suitable for use until the expiration date on the carton when stored and handled as specified. Do not use beyond the expiration date.
- Do not freeze unopened reagent packs.
- Load reagent packs directly from refrigerated storage to minimize condensation.
- Store opened refrigerated reagent packs in a sealed reagent pack storage box that contains dry desiccant.

Calibrator Contents

- 3 sets of VITROS Progesterone Calibrators 1, 2 and 3 (freeze-dried progesterone in human plasma with antimicrobial agent, reconstitution volume 1.0 mL); nominal values 0.0; 4.25 and 120 nmol/L (0.0; 1.33 and 37.7 ng/mL)
- Lot calibration card
- Protocol card
- 24 calibrator bar code labels (8 for each calibrator)

Calibrator Handling

- Use only with reagent packs of the same lot number. Mix thoroughly by inversion and bring to 15–30 °C (59–86 °F) before use. Each pack contains sufficient for a minimum of 6 determinations of each calibrator.
- Handle calibrators in stoppered containers to avoid contamination and evaporation. To avoid evaporation, limit the amount of time calibrators are on the system. Refer to the operating instructions for your system. Return to 2–8 °C (36–46 °F) as soon as possible after use, or load only sufficient for a single determination.

Calibrator Storage and Preparation

Calibrator	Storage Condition		Stability
Unopened	Refrigerated	2–8 °C (36–46 °F)	expiration date
Opened-reconstituted	Refrigerated	2–8 °C (36–46 °F)	≤13 weeks
Opened-reconstituted	Frozen	≤-20 °C (≤-4 °F)	≤13 weeks

- VITROS Progesterone Calibrators are supplied freeze-dried.
- VITROS Progesterone Calibrators are suitable for use until the expiration date on the carton when stored and handled as specified. Do not use beyond the expiration date.
- Reconstitute with 1.0 mL distilled water.
- Opened, reconstituted calibrators may be stored frozen (with no more than 1 freeze-thaw cycle).
- The VITROS Progesterone test uses 25 µL of calibrator for each determination. Transfer an aliquot of each calibrator into a sample container (taking account of the minimum fill volume of the container), which may be bar coded with the labels provided. For details on minimum fill volume of sample cups or containers, refer to the operating instructions for your system.

Specimen Collection, Preparation and Storage

Patient Preparation

No special patient preparation is necessary.

Specimens Recommended

- Serum
- EDTA plasma
- Heparin plasma

Specimens Not Recommended

Do not use turbid specimens. Turbidity in specimens may affect test results.

Special Precautions

IMPORTANT:

Certain collection devices have been reported to affect other analytes and tests.¹¹ Owing to the variety of specimen collection devices available, Ortho Clinical Diagnostics is unable to provide a definitive statement on the performance of its products with these devices. Confirm that your collection devices are compatible with this test.

Specimen Collection and Preparation

- Collect specimens using standard procedures.^{12, 13}
- Thoroughly mix samples by inversion and bring to 15–30 °C (59–86 °F) before use.
- The VITROS Progesterone test uses 25 µL of sample for each determination. This does not take account of the minimum fill volume of the chosen sample container. For details on minimum fill volume of sample cups or containers, refer to the operating instructions for your system.

Handling and Storage Conditions

- Handle samples in stoppered containers to avoid contamination and evaporation.
- The amount of time samples are on the system prior to analysis should be limited to avoid evaporation. Refer to the operating instructions for your system.
- Return to 2–8 °C (36–46 °F) as soon as possible after use, or load sufficient volume for a single determination.
- Serum and plasma samples may be stored for up to 7 days at 2–8 °C (36–46 °F) or 4 weeks at -20 °C (-4 °F).
- Avoid repeated freeze-thaw cycles.

Testing Procedure

Materials Provided

- VITROS Immunodiagnostic Products Progesterone Reagent Pack
- VITROS Immunodiagnostic Products Progesterone Calibrators

Materials Required but Not Provided

- VITROS Immunodiagnostic Products Signal Reagent
- VITROS Immunodiagnostic Products Universal Wash Reagent
- VITROS Immunodiagnostic Products High Sample Diluent A Reagent Pack
- Quality control materials such as VITROS Immunodiagnostic Products RE Controls
- VITROS Immunodiagnostic Products Reagent Pack Storage Box (optional) with desiccant
- Calibrated pipette, distilled water and sample containers for reconstitution of VITROS Progesterone Calibrators

Operating Instructions

Check the inventory regularly to aid the management of reagents and ensure that sufficient VITROS Signal Reagent, VITROS Universal Wash Reagent and calibrated reagent lots are available for the work planned. When performing panels of tests on a single sample, ensure that the sample volume is sufficient for the tests ordered.

For detailed information refer to the operating instructions for your system.

Note:

Do not use visibly damaged product.

Sample Dilution

Serum or plasma (EDTA or heparin) samples with concentrations greater than the measuring range may be automatically diluted on the system up to 10-fold (1 part sample with 9 parts diluent)] by the VITROS Immunodiagnostic and VITROS Integrated Systems with the VITROS High Sample Diluent A Reagent Pack prior to test. Refer to the VITROS High Sample Diluent A Reagent Pack instructions for use.

Default Test Name

The default test name which will appear on patient reports is Progesterone. The default short name that will appear on the test selection menus and laboratory reports is Prog. These defaults may be reconfigured, if required. For detailed information refer to the operating instructions for your system.

Calibration

Calibration Procedure

- Calibration is lot specific; reagent packs and calibrators are linked by lot number. Reagent packs from the same lot may use the same calibration.
- A Master Calibration (a dose response curve covering the full calibration range) is established for each new reagent lot. Concentrations for the linked lot of calibrators are determined from the Master Calibration.
- Ensure that the Master Calibration for each new reagent lot is available on your system.
- Process calibrators in the same manner as samples. Calibration need not be programmed if bar code labels are used; load the calibrators in any order, calibration will be initiated automatically.
- When the calibrators are processed the signal expected for each calibrator is compared against the actual signal obtained. The Master Calibration is then rescaled to reflect the differences between the actual and expected signals. The validity of this calibration curve is assessed against a range of quality parameters, and if acceptable, it is stored for use with any reagent pack of that lot.
- The quality of calibration cannot be completely described by a single parameter. The calibration report should be used in conjunction with acceptable control values to determine the validity of the calibration.
- Recalibration is required after a pre-determined calibration interval, or when a different reagent lot is loaded.
- Calibration results are assessed against a range of quality parameters. Failure to meet any of the defined quality parameter ranges will be coded in the calibration report. For actions to be taken following a failed calibration refer to the operating instructions for your system.

Refer to the operating instructions for your system for detailed instructions on the calibration process.

When to Calibrate

- Calibrate when the reagent pack and calibrator lot changes.
- Calibrate every 28 days.
- After specified service procedures have been performed.
- If quality control results are consistently outside of your acceptable range.

For additional information on when to calibrate, refer to the operating instructions for your system.

Traceability of Calibration

Calibration of the VITROS Progesterone test is traceable to in-house reference calibrators which have been value assigned to correlate to another commercially available test.

Calibration Model

A modified four-parameter logistic curve fit function is used to construct the Master Calibration. The calibration process rescales the Master Calibration to establish a valid stored curve for the VITROS Immunodiagnostic and VITROS Integrated Systems.

Measuring (Reportable) Range

System	Measuring (Reportable) Range
3600 5600 XT 7600 ECi/ECiQ	0.253 [*] –178 nmol/L (0.080–56.0 ng/mL)

^{*} lower limit of measuring range reported by the system software is based on the Limit of Detection.

The lower limit reported by the system can be reconfigured if desired. For details on how to reconfigure the lower limit refer to the operating instructions for your system.

Quality Control

Quality Control Material Selection

VITROS RE Controls are recommended for use with the VITROS Immunodiagnostic and VITROS Integrated Systems. The VITROS RE Controls contain 3 levels of progesterone (low, medium and high). The performance of other commercial control fluids should be evaluated for compatibility with this test before they are used for quality control.

Control materials may show a difference when compared with other progesterone methods if they contain high concentrations of preservatives, stabilizers, or other nonphysiological additives, or otherwise depart from a true human sample matrix.

Appropriate quality control value ranges must be established for all quality control materials used with the VITROS Progesterone test.

Quality Control Procedure Recommendations

- Good laboratory practice requires that controls be processed to verify the performance of the test.
- Choose control levels that check the clinically relevant concentrations.
- To verify system performance, analyze control materials:
 - After calibration
 - According to local regulations or at least once each day that the test is being performed
 - After specified service procedures are performed

If quality control procedures within your laboratory require more frequent use of controls, follow those procedures.

- Analyze quality control materials in the same manner as patient specimens.
- If control results fall outside your acceptable range, investigate the cause before deciding whether to report patient results.
- Refer to published guidelines for general quality control recommendations.¹⁴

For more detailed information, refer to the operating instructions for your system.

Quality Control Material Preparation and Storage

Refer to the manufacturer's product literature for preparation, storage, and stability information.

Results

Results are automatically calculated by the VITROS Immunodiagnostic and VITROS Integrated Systems.

Reporting Units and Unit Conversion

Analyte results are quoted in units of nmol/L or ng/mL. To configure the units, refer to the operating instructions for your system.

Conventional	Alternate
nmol/L (ng/mL × 3.180)	ng/mL (nmol /L × 0.3145)

Limitations of the Procedure

Known Interferences

The VITROS Progesterone test was evaluated for interference consistent with CLSI document EP7.¹⁵ Commonly encountered substances were tested on 3 lots of reagents. The following compounds, when tested, caused the bias shown at the concentration indicated.

Refer to "Specificity" for a list of other compounds tested that did not show interference.

Interferent	Interferent Concentration		Units = nmol/L		Units = ng/mL	
			Analyte Conc.*	Bias**	Analyte Conc.*	Bias**
Dipyrene*	28.5 nmol/L	1000 mg/dL	28.4	5.02	8.93	1.58

* Average concentration of replicate determinations using 2 different lots of reagent.

** Estimate of the average difference observed.

Other Limitations

- The results from this or any other diagnostic test should be used and interpreted only in the context of the overall clinical picture.

INSTRUCTIONS FOR USE

Expected Values and Interpretation of Results

- Heterophilic antibodies in serum or plasma samples may cause interference in immunoassays.¹⁶ These antibodies may be present in blood samples from individuals regularly exposed to animals or who have been treated with animal serum products. Results which are inconsistent with clinical observations indicate the need for additional testing.
- Certain drugs and clinical conditions are known to alter progesterone concentrations *in vivo*. For additional information, refer to one of the published summaries.¹⁷⁻¹⁹
- Biotin levels in serum remain elevated for up to 24 hours after oral or intravenous biotin administration.²⁰

Expected Values and Interpretation of Results

It is recommended that each laboratory establish its own expected values for the population it serves.

Reference Interval

Subject category	Units = nmol/L		Units = ng/mL		n
	Mean	Range	Mean	Range	
Ovulatory Cycle					
Follicular	2.57	0.44–6.47	0.81	0.14–2.03	96
Periovulatory	6.00	1.27–14.2	1.89	0.40–4.47	56
Mid Luteal	39.5	16.6–72.2	12.4	5.22–22.7	66
Luteal	22.1	4.53–52.9	6.95	1.42–16.6	62
Pregnant Females					
1 st Trimester (4–12 weeks gestation)	62.9	20.9–128	19.8	6.57–40.3	137
2 nd Trimester (13–24 weeks gestation)	114	30.7–198	35.9	9.66–62.3	78
3 rd Trimester (25–36 weeks gestation)	288	77.9–1063	90.6	24.5–334	64
Post Menopausal Females	1.27	0.49–3.31	0.40	0.15–1.04	48
Normal Males	2.28	0.66–4.89	0.72	0.21–1.54	50

Progesterone values were determined in samples from a variety of clinical categories. These reference intervals are the central 95% of results in each study.

Interpretation of Results

For patient sample values outside your established reference interval, the system may be configured to display a flag 'LO' or 'HI'. For detailed information refer to the operating instructions for your system.

Performance Characteristics

Limit of Detection

The Limit of Detection (LoD) for VITROS Progesterone is 0.253 nmol/L (0.080 ng/mL), determined consistent with NCCLS document EP17²¹ and with proportions of false positives (α) less than 5% and false negatives (β) less than 1%; based on 698 determinations, with 1 blank and 5 low-level samples. The Limit of Blank (LoB) is 0.085 nmol/L (0.027 ng/mL).

Limit of Blank and Limit of Detection

LoB*		LoD**	
nmol/L	ng/mL	nmol/L	ng/mL
0.085	0.027	0.253	0.080


* Limit of Blank, or the highest value likely to be observed with a sample containing no analyte, replaces the term "analytical sensitivity."

** Proportions of false positives (α) and false negatives (β) were less than 5% and 1% respectively; based on 698 determinations, with 1 blank and 5 low-level samples.

Accuracy (Method Comparison)

Accuracy was evaluated consistent with NCCLS document EP9.²² The plots and table show the results of a method comparison study using patient samples from a variety of clinical categories analyzed on the VITROS ECi/ECiQ Immunodiagnostic System using VITROS Progesterone test lot numbers below 750 compared with those analyzed using the VITROS Progesterone test lot numbers 750 and above. The relationship between the 2 methods was determined by Deming regression.²³

The table also shows the results of method comparison studies²⁴ using patient serum and plasma samples analyzed on the VITROS ECi/ECiQ Immunodiagnostic System compared with those analyzed using the VITROS 3600 Immunodiagnostic System and the VITROS 5600 Integrated System. The relationship between the 2 methods was determined by Deming regression.²³



System	n	Slope	Correlation Coefficient	Conventional Units (nmol/L)		Alternate Units (ng/mL)	
				Range of Samples	Intercept	Range of Samples	Intercept
ECi/ECiQ vs. Comparative Method	210	0.976	0.997	0.370–177	1.23	0.116–55.7	0.387
3600 vs. ECi/ECiQ	106	0.945	0.998	0.91–177	-0.834	0.286–55.7	-0.262
5600* vs. ECi/ECiQ	106	0.949	0.999	0.91–177	-0.410	0.286–55.7	-0.129

* Performance characteristics for the VITROS 5600 System are applicable to the VITROS XT 7600 System.

Precision

VITROS ECi/ECiQ Immunodiagnostic System

Precision was evaluated consistent with NCCLS document EP5.²⁵ Two replicates each of 3 freeze-dried control samples and 1 precision pool were tested on 2 separate occasions per day on at least 20 different days. The experiment was performed using 3 reagent lots on 3 different systems. The data presented are a representation of the product performance.

VITROS 3600 Immunodiagnostic System and VITROS 5600 Integrated System

Precision was evaluated consistent with NCCLS document EP5.²⁶ Two replicates each of 3 freeze-dried control samples were tested on 2 separate occasions per day on at least 20 different days. The experiment was performed using 1 reagent lot on each system. The data presented are a representation of the product performance.

INSTRUCTIONS FOR USE

Prog

Performance Characteristics

System	Units = nmol/L							No. Observ.	No. Days
	Mean Progesterone Conc.	Within-run*		Within-calibration**		Within-lab***			
		SD	CV (%)	SD	CV (%)	SD	CV (%)		
ECi/ECiQ system 1	3.35	0.058	1.8	0.209	6.4	0.297	8.7	88	22
	6.35	0.134	2.1	0.247	3.9	0.300	4.7	88	22
	24.8	0.310	1.2	0.670	2.7	0.750	3.0	88	22
	55.3	0.460	0.8	1.41	2.5	1.61	2.9	88	22
ECi/ECiQ system 2	3.62	0.095	2.5	0.282	7.3	0.363	10.7	88	22
	6.54	0.190	2.8	0.390	5.7	0.472	7.6	88	22
	24.9	0.380	1.5	1.13	4.4	1.26	5.3	88	22
	55.6	0.690	1.2	2.23	3.9	2.23	4.1	88	22
ECi/ECiQ system 3	3.28	0.159	4.9	0.228	7.0	0.241	7.3	88	22
	5.91	0.201	3.4	0.307	5.2	0.350	5.9	88	22
	23.0	0.380	1.7	0.770	3.4	0.800	3.5	88	22
	51.5	0.710	1.4	1.55	3.0	1.59	3.1	88	22
3600	7.82	0.199	2.5	0.320	4.1	0.483	6.1	92	23
	25.0	0.428	1.7	1.09	4.4	1.36	5.4	92	23
	79.9	0.844	1.1	2.63	3.3	2.93	3.7	92	23
5600****	7.10	0.139	2.0	0.234	3.3	0.372	5.2	96	24
	24.3	0.591	2.4	0.850	3.5	1.04	4.2	96	24
	76.7	0.848	1.1	2.16	2.8	2.46	3.2	96	24

* Within-run (repeatability). Between Duplicate precision averaged over all runs

** Within-calibration. Total precision with weighted components of within-run, between-run and between-day variation

*** Within-lab. A measure of the effect of recalibration on total precision, calculated within reagent lot, using data from at least 4 calibrations

**** Performance characteristics for the VITROS 5600 System are applicable to the VITROS XT 7600 System.

System	Units = ng/mL							No. Observ.	No. Days
	Mean Progesterone Conc.	Within-run*		Within-calibration**		Within-lab***			
		SD	CV (%)	SD	CV (%)	SD	CV (%)		
ECi/ECiQ system 1	1.05	0.018	1.8	0.066	6.4	0.093	8.7	88	22
	1.99	0.042	2.1	0.077	3.9	0.094	4.7	88	22
	7.77	0.097	1.2	0.210	2.7	0.235	3.0	88	22
	17.3	0.144	0.8	0.442	2.5	0.505	2.9	88	22
ECi/ECiQ system 2	1.13	0.030	2.5	0.088	7.3	0.114	10.7	88	22
	2.05	0.060	2.8	0.122	5.7	0.148	7.6	88	22
	7.81	0.119	1.5	0.354	4.4	0.395	5.3	88	22
	17.4	0.216	1.2	0.699	3.9	0.699	4.1	88	22
ECi/ECiQ system 3	1.03	0.050	4.9	0.071	7.0	0.076	7.3	88	22
	1.85	0.063	3.4	0.096	5.2	0.110	5.9	88	22
	7.21	0.119	1.7	0.241	3.4	0.251	3.5	88	22
	16.1	0.223	1.4	0.486	3.0	0.498	3.1	88	22
3600	2.46	0.063	2.5	0.101	4.1	0.152	6.1	92	23
	7.86	0.135	1.7	0.343	4.4	0.428	5.4	92	23
	25.1	0.265	1.1	0.827	3.3	0.921	3.7	92	23
5600****	2.23	0.044	2.0	0.074	3.3	0.117	5.2	96	24
	7.64	0.186	2.4	0.267	3.5	0.327	4.2	96	24
	24.1	0.267	1.1	0.679	2.8	0.774	3.2	96	24

* Within-run (repeatability). Between Duplicate precision averaged over all runs

** Within-calibration. Total precision with weighted components of within-run, between-run and between-day variation

*** Within-lab. A measure of the effect of recalibration on total precision, calculated within reagent lot, using data from at least 4 calibrations

**** Performance characteristics for the VITROS 5600 System are applicable to the VITROS XT 7600 System.

Specificity

Substances that do not Interfere

The VITROS Progesterone test was evaluated for interference consistent with CLSI document EP7. ¹⁵ Of the compounds tested, none was found to cause a bias of >10% with the test at the concentrations indicated at progesterone concentrations of 5.03–28.0 nmol/L (1.58–8.80 ng/mL).

Compound	Concentration Tested	
	Bilirubin	0.342 mmol/L
Biotin	81.9 nmol/L	2.00 µg/dL
Hemoglobin	0.31 mmol/L	500 mg/dL
Intralipid	600 mg/dL	6 g/L
Triolein	5.65 mmol/L	500 mg/dL

Cross-Reactivity

The cross-reactivity of the VITROS Progesterone test was evaluated by adding the following substances to control samples containing no progesterone.

Cross-Reactant Tested	Concentration Tested		Mean Value of Cross-reactant Pool		% Cross-reactivity
	µmol/L	ng/mL	nmol/L	ng/mL	
11-Deoxycorticosterone	3.03	1000	126	39.6	4.1
Corticosterone	2.89	1000	95.2	29.9	3.3
17-α-Hydroxyprogesterone	3.03	1000	73.3	23.1	2.4
Testosterone	3.47	1000	5.5	1.73	0.16
Cortisol	6.90	2500	12.4	3.9	0.18
Danazol	29.6	10000	ND*	ND*	ND*
Pregnenolone	3.16	1000	4.8	1.51	0.15
11-Deoxycortisol	2.89	1000	12.3	3.87	0.43
20-α-Dihydroprogesterone	3.16	1000	0.4	0.13	0.01
5-α-Pregnane-3,20-dione	0.316	100	16.5	5.19	5.2
Prednisolone	2.77	1000	ND*	ND*	ND*
Estradiol	3.67	1000	ND*	ND*	ND*
DHEA-SO ₄	25.6	10000	3.7	1.16	0.01
Estriol	3.47	1000	ND*	ND*	ND*
Clomiphene	0.167	100	ND*	ND*	ND*
Bromocryptine	0.133	100	ND*	ND*	ND*

* ND = Not Detectable. Concentration was below 0.25 nmol/L (0.08 ng/mL).

Cross-reactivity was expressed as the mean result obtained for the cross-reactant pool divided by the cross-reactant concentration in percentage term.

$$\% \text{ Cross-reactivity} = \frac{\text{Mean Result for the Cross-reactant Pool}}{\text{Concentration of Cross-reactant}} \times 100$$




































References

1. Siiteri PK, Febres F. In DeGroot et al (Eds). Ovarian Hormone Synthesis - Circulation and Mechanism of Action. *Endocrinology*. 3:1401–1417; 1980.
2. Laufer N et al. The Pattern of Luteal Phase Progesterone and Estradiol in Fertile Cycles. *Am J Obstet Gynecol*. 143:808–813; 1982.
3. Abdulla U et al. Plasma Progesterone levels as an Index of Ovulation. *Br J Obstet Gynaecol*. 90:543–548; 1983.
4. Dunn JF et al. Transport of Steroid Hormones: Binding of 21 Endogenous Steroids to Both Testosterone Binding Globulin and Corticosteroid-binding Globulin in Human Plasma. *J Clin Endocrinol*. 53:58–68; 1981.

5. Hull MGR et al. The Nature of a Single Serum Progesterone Measurement in the Midluteal Phase, as a Criterion of Potentially Fertile Cycle ("Ovulation") Derived from Treated and Untreated Conception Cycles. *Fertil Steril.* 37:355–360; 1982.
6. Hensleigh PA, Fainstat T. Corpus Lutum Dysfunction: Serum Progesterone Levels in Diagnosis and Assessment of Therapy for Recurrent and Threatened Abortion. *Fertil Steril.* 32:396–400; 1979.
7. Radwanska E et al. Plasma Progesterone Levels in Normal and Abnormal Early Human Pregnancy. *Fertil Steril.* 30:398–402; 1978.
8. Summers M et al. Luminogenic Reagent Using 3-Chloro 4-Hydroxy Acetanilide to Enhance Peroxidase/Luminol Chemiluminescence. *Clin Chem.* 41:S73; 1995.
9. CLSI. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline - Fourth Edition.* CLSI document M29-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
10. Regulation (EC) No 1272/2008 of the European Parliament and of the Council of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.
11. Calam RR. Specimen Processing Separator Gels: An Update. *J Clin Immunoassay.* 11:86–90; 1988.
12. CLSI. *Collection of Diagnostic Venous Blood Specimens. 7th ed.* CLSI standard GP41. Wayne, PA: Clinical and Laboratory Standards Institute; 2017.
13. NCCLS. *Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard – Fifth Edition.* NCCLS document H4-A5 [ISBN 1-56238-538-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
14. CLSI. *Statistical Quality Control for Quantitative Measurements: Principles and Definitions; Approved Guideline - Third Edition.* CLSI document C24-A3 [ISBN 1-56238-613-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2006.
15. NCCLS. *Interference Testing in Clinical Chemistry; Proposed Guideline.* NCCLS document EP7-P (ISBN 1-56238-020-6). NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087; 1986.
16. Levinson SS. The Nature of Heterophilic Antibodies and Their Role in Immunoassay Interference. *J Clin Immunoassay.* 15:108–115; 1992.
17. Young DS. *Effects of Drugs on Clinical Laboratory Tests* ed. 4. Washington, D.C.: AACC Press; 1995.
18. Friedman RB, Young DS. *Effects of Disease on Clinical Laboratory Tests.* ed. 3. Washington, D.C.: AACC Press; 1997.
19. Tryding N, Tufvesson C, Sonntag O (eds). *Drug Effects in Clinical Chemistry.* ed. 7. Stockholm: The National Corporation of Swedish Pharmacies, Pharmasoft AB, Swedish Society for Clinical Chemistry; 1996.
20. Scientific Committee on Food. *Opinion of the Scientific Committee on Food on the Tolerable Upper Intake Level of Biotin.* European Commission, SCF/CS/NUT/UPPLEV/55 Final, Brussels, 2001.
21. NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline.* NCCLS document EP17-A (ISBN 1-56238-551-8). CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087, 2004.
22. NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline.* NCCLS document EP9-A (ISBN 1-56238-283-7). NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087, 1995.
23. Deming WE, *Statistical Adjustment of Data.* New York, NY: John Wiley and Sons; 1943.
24. NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline - Second Edition.* NCCLS document EP9-A2 (ISBN 1-56238-472-4). CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2002.
25. NCCLS. *Evaluation of Precision Performance of Clinical Chemistry Devices – Approved Guideline.* NCCLS document EP5-A (ISBN 1-56238-368-X). NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087; 1999.
26. NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition.* NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

Glossary of Symbols

The following symbols may have been used in the labeling of this product.

	Do Not Reuse		Upper Limit of Temperature		Range
	Use by or Expiration Date (Year-Month-Day)		Lower Limit of Temperature		Range of Means
	Batch Code or Lot Number		Temperature Limitation		Midpoint
	Serial Number		Consult Instructions for Use		Revised
	Catalog Number or Product Code		Attention: The Instructions for Use (IFU) has been updated		Supersedes
	Caution		For use in Slide Supply 1		Contains Sufficient for "n" Tests
	Keep Dry (Protect from Moisture/Humidity)		For use in Slide Supply 2		<i>in vitro</i> Diagnostic Medical Device
	Manufacturer		SI Units		Der Grüne Punkt (the Green Dot). Manufacturer follows certain packaging material waste disposal management regulations
	Date of Manufacture		Conventional Units		Estimated within-lab SD
	Authorized Representative in the European Community		Value		Serious Health Hazards
	Corrosive		Flammable		Environmental or Aquatic Toxicity
	Health Hazards		Acute Toxicity		

Revision History

Date of Revision	Version	Description of Technical Changes*
2019-09-06	10.0	<ul style="list-style-type: none"> • Calibrator Contents: updated • Glossary of Symbols: updated • Added EC Representative address

* The change bars indicate the position of a technical amendment to the text with respect to the previous version of the document.

When this Instructions For Use is replaced, sign and date below and retain as specified by local regulations or laboratory policies, as appropriate.

Signature

Obsolete Date

Conditions of supply: all supplies are made subject to the standard terms and conditions of Ortho Clinical Diagnostics or its distributors. Copies of these are available on request.

Distributed in the US by:
Ortho-Clinical Diagnostics, Inc.
100 Indigo Creek Drive
Rochester, NY 14626



Ortho-Clinical Diagnostics
1500 Boulevard Sébastien Brant
B.P. 30335
67411 Illkirch
CEDEX, France



Ortho-Clinical Diagnostics
Felindre Meadows
Pencoed
Bridgend
CF35 5PZ
United Kingdom

VITROS is a trademark of Ortho Clinical Diagnostics.
© Ortho Clinical Diagnostics, 2009–2019

Ortho Clinical Diagnostics