

NABL Accreditation	Specific test	Sample requirement	Instruction for requester	Receipt-Reporting (Turn around time)	Additional examination are accepted upto 24 hrs from the time of collection
Yes	URINE ANALYSIS	Sterile Disposable plastic container; in case of urobilinogen test, it should be in covered container- 10-30 ml Freshly voided preferably morning Urine	Transport at room temperature	6 hours	Not Applicable
No	STOOL EXAMINATION	Sterile Disposable plastic container- Fresh stool	Transport at room temperature	2 hours	Not Applicable
Yes	BODY FLUID EXAMINATION	Sterile disposable EDTA vacuette	Transport at room temperature	6 hours, CSF- 1 hour	Not Applicable for CSF
No	SEMEN ANALYSIS	Sterile Disposable plastic container- Fresh semen	Transport at room temperature, Prior to collection three days of abstinence is advisable	2 hours	Not Applicable